



Application Note

TENSILE, COMPRESSION, IMPACT, MELT FLOW, BRINELL, SHEAR, STRETCH RECOVERY, PUNCTURE, HEAT DISTORTION (VICAT), BEND/FLEXURE, SEAM STRENGTH, TORSION

Proposed Attachment for Skin Testing within Fluid bath



Once the test specimen is mounted on the 'specimen holder' and pinned in place. The main Pin (Pin 1) is extracted to allow the specimen holder to become two parts so that the only component being tested in tension is the skin specimen itself.

Headquarters

Horsham, PA 19044, USA

Tel: +1 215 675 7100

Fax: +1 215 441 0899

Email: info@tiniosolsen.com

web: <http://www.tiniosolsen.com>

International Sales & Service

Salfords, Surrey, RH1 6TL, England

Tel: +44 1737 765001

Fax: +44 1737 765001