

Tinius Olsen Pneumatically Operated Yarn / Fibre Grip HT57



The Tinius Olsen model HT57 is ideal for testing a variety of materials. Examples include;

- **Cords**
- **Filaments**
- **Fibres**
- **Fine Wire**
- **Tyre Cord**
- **Hose Reinforcement.**

Due to the generous radius and large clamping area, invalid tests due to jaw breaks are virtually eliminated.

<i>Maximum capacity</i>	1.5KN (150kgf)
<i>Maximum diameter</i>	4mm
<i>Minimum specimen length at zero jaw separation</i>	150mm
<i>Weight Each</i>	1.4kg
<i>Length Each</i>	150mm

