



MODEL TO-412

Mortar Mixer



This mixer is designed to mix mortars and cement paste to standard requirements and can be operated in either manual or automatic mode.

The mixer features microprocessor control of the speed and mixing program and employs an elliptical mixing motion for thorough and efficient mixing.

TO-412 SPECIFICATIONS

Speed (rpm)	Paddle	Mixing Head
Low	140 ± 5	62 ± 5
High	285 ± 10	125 ± 10
Rated power	180W	
Bowl capacity	5 liters	
Weight	55kg	
Dimension (LxWxH)	760 x 520 x 900mm	

APPLICABLE STANDARDS

- BS 3892-1, 3892-3, 6463-103, 4551-1, ISO 679, EN 196-1, 196-3, 413-3, 459-2, 1744-1, 13279-2, 1015-2, 13395-1, 13454-2

STANDARD FEATURES

- Mixing bowl
- Paddle
- Scraper
- Sand dispenser
- Water burette

ORDERING INFORMATION

- **TO-412-01** Mortar Mixer with sand and water dispenser, 110VAC, 60Hz
- **TO-412-02** Mortar Mixer with sand and water dispenser, 220VAC, 60Hz
- **TO-412-03** Mortar Mixer with sand and water dispenser, 220VAC, 50Hz