


# MODEL TO-418

## Vibration Machine



This custom Vibration Machine is used for vibrating molds with mortar mix at a frequency of  $200 \pm 7$  Hz. The simple design of the machine allows easy assembly and dismantling of the cube molds after vibration. Each machine is certified for its frequency and is supplied with one TO 414 cube mold.



### APPLICABLE STANDARDS

- BS 4550

### ORDERING INFORMATION

- **TO-418-01** Vibration machine with built-in digital timer, 110VAC, 60Hz
- **TO-418-02** Vibration machine with built-in digital timer, 220VAC, 60Hz
- **TO-418-03** Vibration machine with built-in digital timer, 220VAC, 50Hz

### STANDARD ACCESSORIES

- **TO-414** Cube mold 70.6mm

### OPTIONAL ACCESSORIES

- **TO-414** Cube mold 70.6mm
- **TO-417** Cube mold 50mm
- **TO-417-3-CI** 3 gang mold 50mm made from cast iron
- **TO-414-3-NB** 3 gang mold 2inch made from naval brass